

## **AUTOCLAVE COOKER**

### **Aluminium**



## **AUTOCLAVE COOKER**

This vertical laboratory autoclave is designed for safe and efficient steam sterilization of laboratory instruments, glassware, media, and surgical materials. It uses high-pressure saturated steam to ensure complete elimination of bacteria, spores, and microorganisms. The unit features a durable stainless-steel chamber, secure lid-locking system, and a pressure gauge for precise monitoring. Easy to operate and highly reliable, it is ideal for use in research laboratories, hospitals, pharmaceutical units, and educational institutes.

### **FEATURES**

- All Body Made of Aluminium.
- Fitted with pressure gauge.
- The unit has spring loaded safety valve of pressure adjustment and a dead weight safety valve.
- It is hydraulically tested up to 30 PSI & is supplied with Electric fitting, Cord & Plug, Supplied complete with Tripod Stand.
- To work on 220/330 Volt AC Supply.



## TECHNICAL PARAMETERS

| MODEL          | Size (mm) | Size (inches) | Capacity (liter) |
|----------------|-----------|---------------|------------------|
| <b>BRSS-13</b> | 225x275   | 09x11         | 13ltr            |
| <b>BRSS-20</b> | 300x300   | 12x12         | 20ltr            |
| <b>BRSS-25</b> | 300x375   | 12x15         | 25ltr            |
| <b>BRSS-40</b> | 300x500   | 12x20         | 40ltr            |

**BR BIOCHEM LIFE SCIENCES PVT. LTD.**

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: [info@brbls.com](mailto:info@brbls.com) | [www.brbls.com](http://www.brbls.com)



We are available on  **GeM**  
Government  
e Marketplace