

BIOSHAKE



BI-OSI-250 BIOSHAKE

Orbital Shaker-Incubator is a compact bench-top Shaker-Incubator used for mixing of biological liquids as well as for incubation and cultivation of biological liquids according to the operator set program. Built-in microprocessor thermo controller provides constant temperature control in the incubator chamber. It is extremely easy to operate, with very straightforward setup of temperature, speed and time.

FEATURES:

- Instantly display operation information (time, temperature and speed) on the LCD screen.
- Gentle and reliable shaking with direct brushless DC motor
- Running safely -- heating and shaking will be automatically stopped once the lid is open.
- Over-heating protection device is included
- Built-in temperature calibration function.



TECHNICAL PARAMETERS

MODEL	BI-OSI-250
Speed range	40-250rpm
Amplitude	10mm (rotary)
Tray size	Universal tray (Spring clip) 248*248*68(mm)
Input power	220V 50-60HZ
Power	300W
Temperature control range	Room temperature +5°C to 60°C
Temperature setting range	5°C-60°C
Temperature control accuracy	0.1°C
Temperature stability	Plus or minus 0.5C
Timing range	1min-99h59min
Fit vessel height	<210mm
Size Shape	460*325*360mm(LxWxH)
Net weight	18.2kg

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on  **GeM**
Government
e Marketplace