

SAMPLE HOMOGENIZER



BI-SH-6 SAMPLE HOMOGENIZER

Sample Homogenizer is a small size lab instrument, suit for various of samples grinding. Touch screen display, easy to operate: High speed, brings high efficiency. Commonly used for sample pretreatments in industries of biology, chemistry, pharmacology, minerals, medicine, etc.

FEATURES

- Enhanced motor, high quality bearing, touch screen
- High efficiency grinding, stable working, low noise
- Extensive speed range, custom and preset programs
- Transparent lid with protection, easy to use
- Double-position operation

TECHNICAL SPECIFICATIONS

MODEL	BI-SH-6
Speed Range	2450rpm-4450rpm
Sample Capacity	6*2ml /2*5ml
Feed Size	According to tube size
Final Fineness	~5µm
Grinding Beads Diameter	0.1-5mm
Acceleration Time	2s
Deceleration Time	2s
Grinding Method	Dry, Wet
Adapter Material	POM/Nylon
Safety Security	Immediate stop once lid open
Voltage	AC200 ~ 240V 50/60Hz
Dimension(W×D×H)	351mmX215mmX202mm

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, Toll free : 1800 833 6126

EMAIL: info@brbls.com | www.brbls.com



We are available on

