

BACTI-CINERATOR STERILIZER

BR-825A BACTI-CINERATOR STERILIZER

Bacti-cinerator sterilizer is used for incineration and disinfection of small tools and bottles such as inoculation rings, inoculation needles, lab tweezers and triangle and tube mouth. It is easy to operate and widely used in the environment of biological safety cabinets, purification worktables, suction fans and mobile vehicles, etc.

FEATURES:

- Replace alcohol lamp, safe and convenient
- Center temperature up to 825°C, rapid sterilization
- Small and light, easy to use
- Suitable for anaerobic chamber
- Heating chamber material - High-grade ceramic insulation
- Ash organic material to prevent infectious splash and cross contamination
- Automatic shutdown after 3h, extending service life





TECHNICAL PARAMETERS

MODEL	BR825-A
Central High Temperature	825°C+50°C
Standby Temperature	480°C
Heater Head Angle	21°
Max. Sample Diameter	14mm
Heater Head Length	150mm
Power	200W
Weight	1.3kg
Dimension (W×D×H)	105×152×200(mm)
Body	Stainless steel (SS 304)
Heating element	Nichrome coil

BR BIOCHEM LIFE SCIENCES PVT.LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on  **GeM**
Government
e Marketplace