

CHEMIDOC IMAGING SYSTEM



BR-CIS-6 CHEMIDOC IMAGERS

The cute series Chemiluminescence Imagers are the second generation ChemiDoc imagers. The series consist of two models. Both of them adopt super fast cooled digital cameras. Such series are Chemiluminescence detection oriented only.

FEATURES:

- Cute, smart and compact
- Intuitive Software for unlimited users.
- Cooled digital Camera
- Software for unlimited users



TECHNICAL PARAMETERS

MODEL		BR-CIS-6
Standard Parameters	Sensor	Back illuminated SCMOS Sensor
	Resolution	3000x3000
	Quantum Efficiency	Max. 90%
	Temperature	-30 °C regulated
	Dynamic Range	4.80D

Application: Chemiluminescence detection only

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on  **GeM**
Government
e Marketplace