

MICROPLATE SHAKER



BR-MS-100 MICROPLATE SHAKER

This shaker is designed for plates like Elisa plate(96/384 microplate) , cell culture plate (24/48/96 plate), providing adjustable speed for a broad range of applications in immunoassay or dyeing experiments. Safe for use in cool rooms and incubators.

FEATURES:

- LCD display, easy to use
- Equipped with high quality power adapter
- Unique design with four rubbers buffer at corners providing stable shake
- Controlled by PID, precise speed, stable and quiet



TECHNICAL PARAMETERS

MODEL	BR-MS-100
Speed Range	200~1500rpm
Timer	0~99h59m
Amplitude	3mm(Horizontal rotation)
Max. Power	45W
Dimension (W×D×H)	280×278×126 (mm)
Weight	5.6kg
Voltage	AC100~230V 50-60Hz
Power	2.5kw/AC220V.50Hz

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on  **GeM**
Government
e Marketplace