

PCR PLATE CENTRIFUGE



PCR PLATE CENTRIFUGE

BR-PCR-2200 Centrifuge for Microplates is designed for quick spins of samples in PCR plates. Primarily designed for use with 96 or 384 well plates and small volumes, it accepts skirted, non skirted and semi-skirted styles

FEATURES

- Easy to use, real-time display of all operating information.
- Countdown mode support.
- DC brushless motor, low noise, small interference, maintenance-free.
- When opening the cap will automatically stop running.
- Fast acceleration from standstill to maximum speed fully in just 10sec.
- Braking time is short, from the highest speed to a complete standstill in just 4sec



TECHNICAL PARAMETERS

MODEL	BR-PCR-2200
Rotating speed	2200~2800rpm
Timing range	15sec ~1min59sec
Max(RCF)	550g
Acceleration time	≤10sec(0 to 2800rpm)
Brake time	<4sec(2800rpm to 0)
Lid Open Auto Stop	Yes
Short Run	Yes
Capacity	2 PCR plates
Power	45W
Voltage	AC110~220V/50-60Hz
Fuse	250V/3A Φ5×20
Dimension	260x240x226mm
Weight	3.3Kg

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on  **GeM**
Government e Marketplace