

ORBITAL SHAKING INCUBATOR

BRI-111 ORBITAL SHAKING INCUBATOR

An orbital shaker is ideal for a variety of general-purpose shaking applications in cell culture, bacterial growth and suspension, staining and washing procedures.

FEATURES

- The Outer made of M.S. Sheet powder coated & inner chamber is made of thick S.S. Sheet.
- The gap between two walls is fitted with high grade glass wool.
- The glass window makes it possible to inspect specimen without disturbing inner Temperature.
- The unit has a temp. Range from ambient to 60°.
- Temperature is controlled by Dual Display Microprocessor PID Digital
- Temperature Indicator cum Controller with a accuracy of $\pm 0.5^{\circ}$ C.
- Air circulation fan is provided to maintain stability throughout.
- Speed of shaker is variable from 30 RPM to 250 RPM.
- The tray of the shaker is also made of S.S. Sheet with Lotus type clamps to hold flasks.
- Unit is fitted with Digital RPM meter, 999 min digital time.





TECHNICAL PARAMETERS

MODEL	BRI-111		
Capacity	9 flasks of 250 ml	16 flasks of 250 ml	16 flasks of 500 ml
Temp. Range	Ambient to 60°		
RPM Range	Upto 250 RPM		

OPTIONAL (EXTRA COST)

- Timer 0-24 hours.
- 16 flasks of 100 ml instead of 9 flasks of 250 ml.
- 9 flasks of 500 ml or 25 flasks of 100 ml instead of 16 flasks of 250 ml.
- 25 flasks of 250 ml or 36 flasks of 100 ml instead of 16 flasks of 500 ml.

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, Toll free : 1800 833 6126

EMAIL: info@brbls.com | www.brbls.com



We are available on

