

VERTICAL DEEP FREEZER

BRS-350 VERTICAL DEEP FREEZER

A Vertical Deep Freezer, specifically the BRS-350 model, is a 350-liter upright freezer designed for laboratory and storage applications. It features a stainless steel body, six shelves, and a temperature range of -17°C to -24°C . The freezer operates on 240V and 60Hz, making it suitable for use with standard electrical power.

FEATURES:

- Optimal energy saving machines due to better condensing capacity.
- Easy to adjust temperature.
- Uniform Digital temperature display and controller.
- Environment friendly refrigerant and foaming agent.
- Foaming door with dual insulation for better refrigeration effect.
- 7 Numbers of Removable Drawers for easy cleaning and maintenance.



TECHNICAL PARAMETERS

Model No	BRS-350
Capacity(Liters)	350L
Shelves	7 no's shelves
Dimension(WxDxH)	600 x 634 x 1850 (mm)
Temperature range	-17°C to- 24°C
Door Type	Solid Door
Display	Digital Display
Refrigerant	R-134 A
Power	230v/50Hz
Wheels	Adjustable wheels
Automatic Defrost System	Yes
Lockable Door	Yes
No. of Lids /Doors	1

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on

