

CARDBOARD BOX

BR Biochem Cardboard box for 1.5ml 1.8ml Cryogenic Vials Freezer Box are designed for storage of sample as low as -90 centigrade, they can hold both internal and external threaded cryogenic vials as well as storage vials in ultra-low or cyrogenic freezers.

FEATURES

- 81well & 100well available.
- Assorted colors.
- Suitable for use in refrigerators and cryogenic freezer systems temperature down to -80°C .
- Alphanumeric labeling on lid.
- Water proof resistant.
- Compatible with 1.5ml 1.8ml & 2.0ml tubes.
- Made from cardboard.

81well & 100well

PREMIUM QUALITY
★★★



SAFE & RELIABLE

ORDERING INFORMATION

Model	PRODUCT DETAILS	PACK SIZE
BR-CDB-100	BR-CDB-100H Cardboard Cryo Box with hinged lid Pe divider 5/pk	5/PACK
BR-CDB-81	BR-CDB-81H Cardboard Cryo Box with hinged lid PP DIVIDER 5/pk	5/PACK



Premium Quality
Manufactured under stringent quality control standards



Reliable Performance
Consistent results for accurate experimentation



Safe & Durable
Leak-proof, strong and reliable for laboratory use



Customer Support
Dedicated support for all your requirements



1800-833-6126



info@brbls.com



www.brbls.com

