

COOLBOX HOLDER

BR Biochem Coolbox Holder are use with metal coolbox to ensure good cooling capacity, effectively protect the surrounding and bottom of the metal coolbox, prevent slippage and increase stability.

FEATURES

- Silicone material
- Designed for use with Biosharp coolbox
- Protect the bottom of the metal coolbox

METAL COOL BOX

PREMIUM
QUALITY

★★★



SAFE &
RELIABLE

ORDERING INFORMATION

CAT NO.	PRODUCT DETAILS	PACK SIZE
BR032	Metal Cool box with Cool box Holder , Metal Cool box, for 24*1.5ml tubes and 36*0.2ml tubes, Initial storage conditions: -20°C, 4 6h, Provides over 1- 2 hours of 2 to 8°C cooling	1/PK



Premium Quality
Manufactured under
stringent quality
control standards



Reliable Performance
Consistent results
for accurate
experimentation



Safe & Durable
Leak-proof, strong
and reliable for
laboratory use



Customer Support
Dedicated support
for all your
requirements



1800-833-6126



info@brbls.com



www.brbls.com

