

CRYOGENIC GLOVES

BCLS CRYO-S/M/L low temperature of liquid nitrogen protective gloves Material: breathable waterproof composite materials, low-temperature cotton, woven lining, cotton knitted fabrics, Disi.32,38,48,68 cm long optional protection Temperature: -160 ° ~ + 148 ° C



low temperature



FEATURES

- Good heat damage (frostbite, burns) protection, cryogenic liquid splash guard.
- Good multi-layer insulation structure is waterproof, comfortable to wear, light weight, and flexible.
- Good breathability, dry and comfortable inner lining.
- Maintain good comfort when working for a long time in a low temperature environment.
- Machine washable Or use dry cleaning.
- It can be used for liquid oxygen operation in a closed system.

ORDERING INFORMATION

CAT NO.	PRODUCT DETAILS	PACK SIZE
BCLS CRYO-S/M/L	CRYOGENIC GLOVES, Small, Medium & Large	1/Pair



Premium Quality
Manufactured under stringent quality control standards



Reliable Performance
Consistent results for accurate experimentation



Safe & Durable
Leak-proof, strong and reliable for laboratory use



Customer Support
Dedicated support for all your requirements



1800-833-6126



info@brbls.com



www.brbls.com



FOLLOW US

