

ICE FLAKING MACHINE



ICE FLAKING MACHINE

- High-quality stainless steel is used for outer body and they are anti-corrosive and durable, and the all-in-one independent structure makes them simple and compact and space-saving. Bin is made of good quality plastic.
- Their boxes are protected with fluorine-free foam thermal cover featuring good insulation effect. Their internal bladders are of fluorine-free bacterial inhibitory-type. They are both energy-saving and environmental friendly.
- High quality and efficient use of fluorine-free R134a compressors, parts and components have passed relevant safety certifications, and they are safe and reliable.
- They have ice full display, water deficiency display, overcooling protection display and malfunction warning display, etc. for shutdown protection. The stop automatically in case of ice full or water deficiency and start automatically when power or water supply is resumed.
- The ice made is small shaped, amorphous granular flake chipped ice, which can infiltrate into narrow space, featuring fast cooling speed, good ice bathing effect. It is specially designed for laboratories.
- There are power switch and function indicators deployed in the front with detailed instructions. It is intuitive and easy to use. Various safety performance indicators have passed electrical tests, they are safe and reliable.
- Flake shaped ice is widely used in hospitals, schools, laboratories, research institutes and so on, as well as for food preservation in supermarkets, fishing refrigeration, medical applications, chemical, food processing and other industries. It has wide range of applications.



BR BIOCHEM
LIFE SCIENCES PVT. LTD.

TECHNICAL PARAMETERS



Item no	BRIF-20	BRIF-30	BRIF-50	BRIF-70	BRIF-85	BRIF-100	BRIF-120	BRIF-150	BRIF-200	BRIF-250	BRIF-300
Ice Making Capacity (kg/24h)	20	30	50	70	85	100	120	150	200	250	300
Ice Storage Capacity (kg)	10	10	15	25	25	40	40	40	60	60	150
Input Current (w)	280	180	380	460	480	520	550	680	805	1370	1360
Ambient Temperature	10~25℃										
Water Temperature	10℃										
PowerSupply	AC220V±10%, 50/60HZ; 110V±10%,60HZ										
Refrigeration Type	Air cooling										
Tank shell	Plastic										
Compressor/ Refrigerant	CFC-free/R134a										
Standard Accessories Water	Water inlet pipe*1, Drainage pipe*1, Ice spoon*1										
Optional Accessories	Filter (Water pressure needs to be within 1.5~3Kg/cm2)										
Net Weight (kg)	28	32	36	57	57	58	58	60	85	91	96
Gross Weight (kg)	36	36	39	62	62	63	63	65	91	98	105
External Size (WxDxH) mm	330 x 470 x 605		400 x 543 x 720		500 x 611 x 940		500 x 611 x 940		560 x 800 x 1750		
Package Size (W x D x H) mm	405 x 555 x 665		470 x 615 x 780		570 x 685 x 1000		570 x 685 x 1000		650 x 850 x 1750		
Applicable Area	Various types of experiments, supermarkets, hotels and other seafood products stores etc.										

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on

