

MICRO-CENTRIFUGE TUBES



BULK



BR Biochem Micro-centrifuge Tubes are made of transparent polypropylene material. Micro-centrifuge tubes are available in a variety of volumes. Mainly used for sample storage and centrifugation. All tubes feature frosted writing surfaces on the cap and side for smear resistant labeling.

FEATURES

- Manufactured from polypropylene material
- Positive flat cap seal design
- With graduation and marking area
- Withstand temperature down to -80°C
- Autoclavable
- Bulk pack

ORDERING INFORMATION

CAT NO.	PRODUCT DETAILS	PACK SIZE
B0.5MCT	Centrifuge Tube 0.5ml, Bulk Pack	1000/PK
B1.5MCT	Centrifuge Tube 1.5ml, shield lid, Bulk Pack	1000/PK
B2MCT	Centrifuge Tube 2ml, Bulk Pack	1000/PK



Premium Quality
Manufactured under
stringent quality
control standards



Reliable Performance
Consistent results
for accurate
experimentation



Safe & Durable
Leak-proof, strong
and reliable for
laboratory use



Customer Support
Dedicated support
for all your
requirements



1800-833-6126



info@brbls.com



www.brbls.com



FOLLOW US

