

MINI CENTRIFUGE



MINI 10K PLUS MINI CENTRIFUGE

The mini centrifuge is novel in appearance and versatile. It is equipped with two centrifugal rotors and multiple test tube sleeves for 2.0ml/1.5ml/0.5ml/0.2ml centrifuge tubes and 8 × 2 × 0.2ml PCR tube strips. Fully transparent round top cover with multi rotor equipped, and an open cover pause function and electronic timing fully embodied the design concept of humanization.

FEATURES:

- Large radius with 8-hole rotor design, 1.5-fold stronger centrifugal force
- Top cover opens automatically for easy one-handed operation
- Snap-on design for easy rotor replacement and removal
- Noise reduction function, near silent operation
- The maintenance-free motor is durable, safe and reliable
- Exquisite appearance and compact size



TECHNICAL PARAMETERS

MODEL	MINI 10 K PLUS
Speed	3000-10000rpm
Max. RCF	500-5400g
Capacity	Centrifuge Tube: 8×2.0ml/1.5ml/0.5ml/0.2ml PCR tube strips: 8x2x0.2ml
Timer	1~99h59m
Noise Level	≤55dB
Dimensions (WxDxH)	160×176×135 (mm)
Weight	≤1.5kg
Voltage	220V/110V 50-60Hz

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on

