

VERTICAL ELECTROPHORESIS

MP-8003 VERTICAL ELECTROPHORESIS

Vertical electrophoresis tank is suitable for biological research, the separation, purification, preparation of nucleic acid and protein samples, etc., the analytical protein electrophoresis can meet the purity identification, the analysis of complex protein samples, and also suitable for nucleic acid electrophoresis. The electrophoresis tank uses a more scientific design and diverse functions. It can flexibly choose 1-4 gels for experiments. Simple operation, stability and safety are very friendly to the experimenter, and can prevent misassembly.

FEATURES:

1. Truly compatible with domestic and foreign brands of prefabricated glue.
2. Product parameters, and accessories are compatible with imported products.
3. The edge cushion is fixed on the long glass plate, which effectively improves the accurate alignment of the glass plate and reduces the glue leakage.
4. The rubber frame of the cam card lock is easy to operate, and it can accurately align the glass plate in any plane.
5. Special plastic Well Combs do not inhibit the gel polymerization reaction, during the built-in ridges, can effectively improve the uniform gel polymerization.
6. High-purity platinum electrode wire, to achieve a better conductive performance;
7. It can accommodate 1-4 pieces of prefabricated glue and complete up to 4 small prefabricated glue electrophoresis within 1 hour.



TECHNICAL PARAMETERS

MODEL	MP-8003
Gel number	1-2 block
Gel thickness	1.0mm, 0.75mm, 1.5mm
No. of samples	10 teeth optional 1.5 teeth;
Gel size	Hand glue (83×73cm)
Prefabricated glue	8.6×6.8cm

ACCESSORIES	QUANTITY
Main Body 4-Gels	1
Comb	5
	5
Spacer Glass Flat	5
Short Glass Flat	10
Gel Base	4
Gel-Casting Frame	4
Hinge	8
Sealing Gasket	5
U-Sealing Strip	4
Side Fastener	4
Electrode Head Electrophoresis Core	1
Mushroom Head Electrophoresis Core	1
Buffer Dam	1
Gel Shovel	5
Sample-Adding Frame	1
Upper Lid (with Cable)	1
Lower Lid	1
Power Cable	1
Packing List	1

BR BIOCHEM LIFE SCIENCES PVT.LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com | www.brbls.com



We are available on 