

PCR TUBES & STRIP


PCR Tubes & 8 Strip Tubes are designed to fit a wide variety of thermal cyclers. Ultra-thin and consistent wall thickness allows precise thermal transfer, ensures a tight fit to prevent evaporation.

FEATURES

- Manufactured from polypropylene material
- Designed for effective thermal conductivity
- Suitable for all kinds of universal PCR reactions
- Flat cap and domes cap single tubes are available
- Autoclavable
- Non-sterile

BULK

PREMIUM
QUALITY
★★★



SAFE &
RELIABLE

ORDERING INFORMATION

CAT NO.	PRODUCT DETAILS	PACK SIZE
BL60-0002	PCR tube 0.2ml, Optical clear with optical flat cap	1000/PK
BL60-008	8-pc Tube Strip, Vol.0.2ml, Tube and Lid, Individual Packing	125/PACK
BL60-008A	8-pc Tube Strip, Vol.0.2ml, Tube attached cap	120/PACK



Premium Quality
Manufactured under
stringent quality
control standards



Reliable Performance
Consistent results
for accurate
experimentation



Safe & Durable
Leak-proof, strong
and reliable for
laboratory use



Customer Support
Dedicated support
for all your
requirements



1800-833-6126



info@brbls.com



www.brbls.com



FOLLOW US

