

ULTRASONIC CLEANER DEGASSING

ULTRASONIC CLEANER DEGASSING

The ultrasonic cleaning machine is made of high quality stainless steel , which has strong corrosion resistance and long service life. Using original Korean high performance and high quality ultrasonic transducer, with advanced bonding technology, high electro-acoustic conversion efficiency and strong ultrasonic output power. With automatic temperature heating device, temperature control range: room temperature ~ 50 °C. Suitable for parts cleaning in various industries, such as precision electronic components, eyeglasses and watch parts, optical clock parts, hardware machinery parts, jewelry, semiconductor silicon wafers, polyester filter cores/spray plates, electronics, electrical appliances, hardware, machinery, communications, etc. Cleaning, cleaning and cleaning of enterprises, institutions, and commercial fields such as biology, medical, optics, jewelry, watches, chemicals, textiles, automobiles, ships, aviation, and mining, etc., cleaning of medical devices, and cleaning of parts before and after electroplating.





TECHNICAL PARAMETERS

MODEL		BR-USC-4L	BR-USC-6L	BR-USC-10L	BR-USC-15L	BR-USC-22L	BR-USC-30L
Oscillation	Power	180 W	180 W	220W	350 W	480 W	600 W
	Frequency	40 kHz	40 kHz	40 kHz	40 kHz	40 kHz	40 kHz
	Timer Range	1-30 min.	1-30 min.	1-30 min.	1-30 min.	1-30 min.	1-30 min.
Heating	Power	100W	200W	260 W	300W	500W	600 W
	Max.Temp.	176°F 80°C	176°F 80°C	176°F 80°C	176°F 80°C	176°F 80°C	176°F 80°C
Tank	Capacity	3.2L	5.8L	9.6L	14.5L	21.5L	29L
	Total Volume	4L	6L	10L	15L	22L	30L
	Dimensions	26.5x15.4x15 cm	30x15.4x15 cm	30x24.3x15 cm	33x30.3x15 cm	50.5x30.3x15 cm	50.5x30.3x20 cm
Overall Dimensions		28x17.7x28.1cm	32.5x17.7x28.1cm	32.7x26.7x27.4cm	35.5x32.5x27cm	53.5x33.3x28 cm	53.5x33.3x32.5cm
Net Weight		3kg	4.6kg	5.8 kg	7.8kg	10.6 kg	12.2kg

BR BIOCHEM LIFE SCIENCES PVT. LTD.

06, S.F. Gopal Nagar , Fateh Nagar, New Delhi-110018 ,India

Ph: +91 11 43033662/4, +91 8860082777

EMAIL: info@brbls.com |

www.brbls.com



We are available on  **GeM**
Government e Marketplace