

Cat. No. (SKU)	Product Name	Size Range (bp)	No. of Bands	Tracking Dyes / Notes	Pack Sizes Available	Storage
BM014-10	10 bp DNA Ladder	80 - 300 bp	23	Rec. Load: 2-5 µl/Lane, Conc: 168 ng/µl	500 µl / 10x500 µl	4°C / - 20°C
BM014-20	20 bp DNA Ladder	60 - 300 bp	13	Rec. Load: 2-5 µl/Lane, Conc: 128 ng/µl	500 µl / 10x500 µl	4°C / - 20°C
BM012-R500	50 bp DNA Ladder	50 - 1,500 bp	17	Orange G tracking dye	500 µl / 10x500 µl	4°C / - 20°C
BM001-R500	100 bp DNA Ladder	100 - 3,000 bp	12	Orange G & Xylene Cyanol FF	500 µl / 10x500 µl	4°C / - 20°C
BM010-R500	1kb bp DNA Ladder	250 - 10,000 bp	13	Orange G & Xylene Cyanol FF	500 µl / 10x500 µl	4°C / - 20°C
BM011-R500	kPlus DNA Ladder	100 - 10,000 bp	19	Orange G & Xylene Cyanol FF	500 µl / 10x500 µl	4°C / - 20°C
BM013-R500	XLarge DNA Ladder	250 - 25,000 bp	14	Wide range ladder	500 µl / 10x500 µl	4°C / - 20°C