

Cat. No. (SKU)	Product Name	Key Application / Description	Pack Sizes	Storage Temp
BCMB203-0100	All-IN-ONE PCR Master Mix	Direct gel loading, no staining required	100 rxn / 10x100 rxn	-20°C
BCQ050 / BCQ500	Taq DNA Polymerase (5U/μl)	Thermostable DNA amplification up to 5 kb	500 U / 5x500 U	-20°C
BCN-1	dNTP Mix (10 mM each, SBS)	Ready-to-use dATP, dCTP, dGTP, dTTP (≥99% HPLC)	1 mL / 10x1 mL	-20°C
BCGB100 / BCGB500	cDNA Synthesis Kit	Sensitive M-MLV RTase solution for RT-PCR	100 rxn / 500 rxn	-20°C
BCLD005	Gel Loading Dye (6X)	Premixed dye with Bromophenol Blue & Xylene Cyanol	5 x 1 mL	Room Temp
HGV-II	GoodView™ Nucleic Acid Stain	Non-carcinogenic EtBr alternative for dsDNA/RNA	1 mL / 5x1 mL	2 to 8°C
BR100	Trizol Reagent	Total RNA purification from tissue, cells, bacteria	100 mL / 5x100 mL	2 to 8°C
BCB1125	T4 DNA Ligase (5U/μl)	Catalyzes phosphodiester bonds for cohesive/blunt DNA	500 U / 2500 U	-20°C